



HAMILTON-G5

The modular high-end ventilation solution



The modular high-end ventilation solution

HAMILTON-G5 - Intelligent Ventilation built in

The HAMILTON-G5 is Hamilton Medical's most modular high-end mechanical ventilator. A large set of standard features and options allows you to tailor the HAMILTON-G5 to your needs. The integrated cuff pressure controller, IntelliCuff®, continuously monitors and maintains a set cuff pressure in real-time. The Protective Ventilation Tool (P/V Tool Pro) provides a repeatable method to assess recruitability, and to perform recruitment maneuvers.

The HAMILTON-G5 also supports transpulmonary pressure measurement.

In addition to many other features, the HAMILTON-G5 includes Hamilton Medical's standardized Ventilation Cockpit user interface and the intelligent ventilation mode, Adaptive Support Ventilation (ASV®).

Your benefits

- ✓ Advanced ventilation modes including ASV and INTELLiVENT-ASV
- ✓ Automated cuff pressure controller IntelliCuff
- ✓ P/V Tool Pro for lung assessment and recruitment
- ✓ Transpulmonary pressure measurement
- ✓ High flow oxygen therapy
- ✓ Adult, pediatric, and neonatal ventilation



“



I think the HAMILTON-G5 is a unique machine, and it is almost revolutionizing the way we, as critical care physicians, take care of our patients.

Dr. Shaikh A. Hai, FACS, Chief of Surgical Critical Care
The Methodist Hospital, Houston (TX), USA

Product overview

- 1 Flexible device configuration
- 2 Patient interfaces and ports
- 3 Press & turn control knob
- 4 Ventilation Cockpit
- 5 360° visible alarm lamp
- 6 Option modules
- 7 Remote control of the humidifier from the Ventilation Cockpit
- 8 Multimedia interfaces







“

What I like the most about the HAMILTON-G5s is the monitoring parameters and the ability to trend those up to 96 hours. I have been able to use that on a disease-specific basis and trend data that I couldn't do before.

Craig Jolly, RRT, Adult Clinical Education Coordinator
University Medical Center, Lubbock (TX), USA



Ease of use

Intuitive operation

In close cooperation with users and ventilation experts, our engineers have designed the HAMILTON-G5's user interface to allow intuitive operation and direct access to important settings. All Hamilton Medical ventilators are operated according to the same principles, which makes switching between different devices very easy.

Easy-to-understand monitoring

Ventilators display large amounts of data that is often difficult to interpret. The configurable touch screen display, referred to as the Ventilation Cockpit, consolidates the diverse monitoring data, and presents it numerically and in various graphic panels. These easy-to-understand views provide an at-a-glance overview of the patient's current ventilation status, and offer a reliable basis for therapy decisions.

More time for your patients

In ASV and INTELLiVENT-ASV modes, the ventilator continuously adjusts to the patient's lung condition and breathing activity. This means fewer user interactions are required¹⁾ and fewer alarms are generated²⁾, giving you more time for your patients.

1 Beijers AJ. Intensive Care Med. 2014 May;40(5):752-3. | 2 Zhu F. Anesthesiology. 2015 Apr;122(4):832-40.

More safety and comfort for your patients

Lung-protective ventilation

Hamilton Medical's intelligent ventilation modes, ASV and INTELLiVENT-ASV, automatically employ lung-protective strategies to minimize complications such as AutoPEEP and volutrauma/barotrauma. Within the rules of this lung-protective strategy, they encourage the patient to breathe spontaneously.

According to several studies, these intelligent ventilation modes help reduce ventilation time in various patient groups. They decrease work of breathing and improve patient-ventilator synchrony.

P/V Tool Pro offers you additional support to carry out your lung-protective ventilation strategy.

ASV decreases ventilation time


Publications show that

- ASV supports earliest possible spontaneous breathing by the patient ^{1), 2)}
- ASV shortens ventilation time in various patient groups ^{1), 2)}

INTELLiVENT-ASV is safe and efficient

Publications show that

- INTELLiVENT-ASV is first among all modes in terms of capabilities related to safety, comfort, and weaning ³⁾
- INTELLiVENT-ASV was safe and efficient for post-cardiac surgery patients, and requires fewer interactions than conventional modes ⁴⁾



Patients in the medical intensive care unit could be extubated earlier following the introduction of ASV.⁵

Time to Extubation
Readiness (d)

Increase in efficiency

Integrated commercial considerations

Ventilators are capital goods that need to be evaluated for cost efficiency. Factors such as treatment costs and the use of human resources play an important role in this process. Assembled with an extensive standard equipment package that is easy to maintain, Hamilton Medical ventilators are an attractive investment with respect to purchase price and operating costs.

Reduction of treatment costs

Each eliminated ventilation day significantly reduces treatment costs – on average by 1,500 USD.¹⁾ It has been shown that the use of Hamilton Medical ventilators and ASV can reduce ventilation time. In addition, the ventilator is now available for the next patient much earlier. A shorter ventilation time also reduces the risk of ventilator associated pneumonia (VAP), which may result in costs of up to 57,000 USD per case.²⁾

Improved use of human resources

Hamilton Medical ventilators, along with ASV, can reduce the time needed for standard settings and alarm management while maintaining ventilation quality.^{3), 4)} This frees up time for other aspects of patient care. Thanks to ease of operation, consistent operating concepts across devices, and the free e-learning offerings from Hamilton Medical, the effort for education and training is also reduced.

One device for all applications

The HAMILTON-G5 makes life easier for caregivers by supporting a wide range of ventilation therapies from invasive and noninvasive to high flow oxygen therapy. The patient can stay on the same ventilator as long as respiratory support is needed.

Intelligent ventilation solutions

ASV and INTELLiVENT-ASV: automated, lung-protective ventilation modes

The HAMILTON-G5 features the intelligent ventilation mode ASV. ASV maintains an operator-set minute volume, and continuously determines respiratory rate, tidal volume, inspiratory pressure, and inspiratory time depending on the patient's lung mechanics and effort. ASV adapts ventilation breath-by-breath, 24 hours a day, from intubation to extubation.

The unique INTELLiVENT-ASV mode takes the proven concepts of ASV a step further: the clinician sets targets for PetCO_2 and SpO_2 for the patient. INTELLiVENT-ASV then automates the ventilator controls to reach these targets while taking into account physiologic inputs from the patient (PetCO_2 , SpO_2 , lung mechanics, spontaneous breaths). It automatically sets controls relevant to oxygenation (PEEP, oxygen) and ventilation (respiratory rate, tidal volume, inspiratory pressure, and inspiratory time).


Both ASV and INTELLiVENT-ASV automatically employ lung-protective strategies to minimize complications from AutoPEEP and volutrauma/barotrauma. Within the rules of this lung-protective strategy, ASV and INTELLiVENT-ASV encourage the patient to breathe spontaneously. In addition, INTELLiVENT-ASV provides a configurable automated weaning protocol (Quick Wean) including spontaneous breathing trials.

Protective ventilation with P/V Tool Pro

P/V Tool Pro allows an easy assessment of lung recruitability, and provides a simple and safe way to perform lung-recruitment maneuvers. In combination with transpulmonary pressure measurement, P/V Tool Pro allows you to more precisely assess lung recruitability and perform recruitment maneuvers in ARDS patients.



Ventilation Cockpit with INTELLiVENT-ASV



P/V Tool Pro window

Neonatal ventilation

Tidal volumes as low as 2 ml

With the neonatal option, the HAMILTON-G5 provides tidal volumes as low as 2 ml for an effective, safe, and lung-protective ventilation of even the smallest preterm infants.¹⁾ The proximal flow sensor, specifically developed for neonates, precisely measures the pressure, volume, and flow directly at the infant's airway opening, ensuring the required trigger sensitivity. This provides improved synchronization and less work of breathing.

Adaptive synchronization, even with uncuffed tubes

Leaks are one of the issues encountered in the ventilation of neonates, as a result of using uncuffed tubes. The IntelliTrig leak compensation function automatically adjusts the inspiratory and expiratory trigger sensitivity to potential leaks. This enables adaptive synchronization with the neonate's breathing pattern.

nCPAP - Automatic adaptation, fewer interventions

The HAMILTON-G5's nCPAP mode is designed in such a way that you only need to set the desired CPAP pressure. The flow is subsequently adjusted automatically based on the patient condition and potential leaks. This prevents unintended peak pressures and guarantees highly efficient leak compensation. Flow adjustment occurs very rapidly due to near-patient pressure measurement and the high sensitivity of the measurement.



Neonatal ventilation with HAMILTON-G5



Fewer interventions and increased safety with nCPAP

1 Volume-targeted versus pressure-limited ventilation in the neonate (Review), 2011 Morley CJ.

Tailor the ventilator to your needs




Automated ventilation with INTELLiVENT-ASV (optional)

INTELLiVENT-ASV is based on the proven Adaptive Support Ventilation mode. It continuously monitors patient conditions, and safely adjusts parameters to keep the patient within target ranges, with minimal clinician interaction, from intubation to extubation. INTELLiVENT-ASV also provides tools to promote early, automated weaning (Quick Wean).




High flow oxygen therapy (optional)

To accommodate your patient's needs, the HAMILTON-G5 allows you to switch easily between any ventilation mode and high flow oxygen therapy. Using the same breathing circuit, you simply need to change the ventilation mode and the patient interface.



Assess lung recruitability and perform lung recruitment


P/V Tool Pro provides a fully automatic maneuver that records the static pressure-volume curve quickly and easily at the bedside. P/V Tool Pro allows an easy assessment of lung recruitability, and provides a simple and safe way to perform lung-recruitment maneuvers.



Transpulmonary pressure measurement

The ventilator provides an auxiliary port to connect an esophageal balloon catheter and displays esophageal (Pes) and transpulmonary pressures (Ptranspulm). Transpulmonary pressure can be used in combination with the P/V Tool Pro for assessing recruitability and performing recruitment maneuvers.

Tailor the ventilator to your needs



Continuous cuff pressure management (optional)

IntelliCuff is a new noninvasive continuous cuff pressure controller integrated with the ventilator. IntelliCuff continuously monitors and adjusts set cuff pressure in real-time.



Volumetric capnography (optional)

Proximal flow and CO₂ measurement enables the HAMILTON-G5 to perform up-to-date volumetric capnography. This provides an important basis for the assessment of ventilation quality and metabolic activity. Alternatively the HAMILTON-G5 also offers sidestream CO₂ monitoring.



Integrated Aerogen® nebulizer (optional)

An integrated synchronized Aerogen nebulizer helps to conserve expensive medications. It improves drug delivery efficiency, and offers the potential to reduce drug and personnel costs associated with in-patient treatment, while maintaining the integrity of ventilator-dependent care.



SpO₂ measurement (optional)

The integration of the SpO₂ option with Hamilton Medical ventilators offers a new and innovative solution for noninvasive measurement. Hamilton Medical also provides a comprehensive portfolio of SpO₂ sensors.

Hamilton Medical

Intelligent Ventilation since 1983

In 1983 Hamilton Medical was founded with a vision: To develop intelligent ventilation solutions that make life easier for patients in critical care and for the people who care for them. Today, Hamilton Medical is a leading manufacturer of critical care ventilation solutions for a wide variety of patient populations, applications, and environments.

The right ventilation solution for any situation

The ventilators from Hamilton Medical ventilate all of your patients; in the intensive care unit, during an MRI procedure and in all transport situations, from the neonate to the adult. Each of these ventilators is equipped with the same standardized user interface and uses the same Intelligent Ventilation technologies. This enables Hamilton Medical ventilators to help you to

- ✓ Increase the comfort and safety of your patients
- ✓ Make life easier for the caregivers
- ✓ Increase efficiency and return on investment

