


Navigate with confidence.



Roadrunner® PC
HYDROPHILIC WIRE GUIDE



Navigate with confidence.

Navigating anatomy can be difficult. Choosing the right wire guide can help. The Roadrunner PC Hydrophilic Wire Guide can be used to gain ureteral access, to establish a tract, and to assist in the placement, replacement, and exchange of medical devices during urological procedures. It may also be used for catheter positioning and exchange in a tortuous or kinked ureter, for traversing a stone en route to the kidney, or in cases demanding enhanced control and high radiopacity.


Features and benefits

- Hydrophilic coating for reduced friction as the wire advances
- Platinum tip for visualisation during fluoroscopy
- Nitinol core for kink resistance during placement
- One-to-one torque for control as you place the wire

Options


- Body configurations: Standard or stiff
- Tips: Angled, straight, or double-flexible
- Hydrophilic coatings: AQ® or Slip-Coat®
- Ink marks: A total of 10 ink marks, one every 5 cm, starting at the distal tip (on select wire guides)

The Roadrunner wire guide is packaged with an Olcott Torque Device.



Roadrunner® PC

HYDROPHILIC WIRE GUIDE



Order Number	Reference Part Number	Diameter inch	Length cm	Tip Configuration	Coating
Standard-Body Wire Guides					
G17540	RFSPC-035145	0.035	145	3 cm angled	Slip-Coat
G17541	RFSPC-038145	0.038	145	3 cm angled	Slip-Coat
G17542	RFSPC-035145-0	0.035	145	3 cm straight	Slip-Coat
G17539	RFSPC-038145-0	0.038	145	3 cm straight	Slip-Coat
Standard-Body Wire Guides with Ink Marks					
G18629	RFSPC-035145-0-I-AQ	0.035	145	3 cm straight	AQ
Stiff-Body Wire Guides*					
G34132	RPC-035145-5	0.035	145	8 cm angled	
G34130	RPC-038145-5	0.038	145	8 cm angled	Slip-Coat
G34131	RPC-035145-0-5	0.035	145	8 cm straight	Slip-Coat
G34129	RPC-038145-0-5	0.038	145	8 cm straight	Slip-Coat
Stiff-Body Wire Guides – Double Flexible tips with Ink Marks					
G18792	RSPC-035145-DF-I-AQ	0.035	145	5 cm double flexible	AQ
G18794	RSPC-038145-DF-I-AQ	0.038	145	5 cm double flexible	AQ

*Sold in boxes of 5

Slip-Coat is a registered trademark of Angiotech Biocoatings Corp.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



1CC, URO-A4