




Cook® EchoTip® Double Lumen Aspiration Needle

Used for aspiration and flushing of oocytes from ovarian follicles.


- The double-lumen design allows simultaneous or intermittent flushing of media and aspiration of follicular fluid.
- EchoTip enhances the visualization of the needle tip under ultrasound.
- The standard aspiration line is 75 cm long.



Order Number	Reference Part Number	Needle gage	Needle Length cm	Flushing Line Length cm
A Bevel				
G30408	K-OPSD-1630-A-S	16	30	70
G55786	K-OPSD-1630-A-S-US	16	30	70
G30409	K-OPSD-1630-A-L	16	30	100
G55785	K-OPSD-1630-A-L-US	16	30	100
G30414	K-OPSD-1635-A-S	16	35	70
G55788	K-OPSD-1635-A-S-US	16	35	70
G30415	K-OPSD-1635-A-L	16	35	100
G55787	K-OPSD-1635-A-L-US	16	35	100
G55795	K-OPSD-1730-A-S-US	17	30	70
G55794	K-OPSD-1730-A-L-US	17	30	100
G55797	K-OPSD-1735-A-S-US	17	35	70
G55796	K-OPSD-1735-A-L-US	17	35	100
B Bevel				
G55791	K-OPSD-1630-B-S-US	16	30	70
G30411	K-OPSD-1630-B-L	16	30	100
G55790	K-OPSD-1630-B-L-US	16	30	100
G55792	K-OPSD-1635-B-S-US	16	35	70
G55798	K-OPSD-1735-B-S-US	17	35	70

Use reference part numbers with the suffix "-US" if ordering in the U.S., Canada, or Latin America.


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Immature Ovum Aspiration Needle Set

Used for laparoscopic or ultrasound-guided transvaginal aspiration of oocytes from small ovarian follicles.

- The EchoTip® design extends around the bevel, enhancing the visualization of the needle tip under ultrasound.
- The guide needle adds rigidity to the aspiration needle and is used to puncture the vaginal wall.
- The stopper design accepts the tip of a Luer slip syringe, facilitating follicle flushing.




Order Number	Reference Part Number	Guide Needle gage	Guide Needle Length cm	Aspiration Needle gage	Aspiration Needle Length cm	Aspiration Line Length cm
G38609	K-IOPS-2035-1730	17	30	20	35	60

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Otrieva® Tapered Ovum Aspiration Needle

Used for laparoscopic or ultrasound-guided transvaginal aspiration and flushing of oocytes from ovarian follicles.



Order Number	Reference Part Number	Needle gage	Length cm
A Bevel			
G49109	K-TIVM-172035	17 tapers to 20	35
G56541	K-TIVM-172035-US	17 tapers to 20	35
G49108	K-TIVM-172030	17 tapers to 20	30
G56542	K-TIVM-172030-US	17 tapers to 20	30

Use reference part numbers with the suffix "-US" if ordering in the U.S., Canada, or Latin America.


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle

Used for laparoscopic or ultrasound-guided transvaginal aspiration and flushing of oocytes from ovarian follicles.

- The ergonomic handle provides comfort and control during use.
- The extra-stiff, smooth needle cannula enables precise placement of the needle tip in the follicle.
- EchoTip enhances the visualization of the needle tip under ultrasound.



Order Number	Reference Part Number	Needle gage	Needle Length cm	Aspiration Line Length cm
A Bevel				
Heavy Wall				
G30431	K-OSN-1635-A-60	16	35	60
G30459	K-OSN-1635-A-90	16	35	90
G30437	K-OSN-1730-A-60	17	30	60
G29620	K-OSN-1735-A-60	17	35	60
G29624	K-OSN-1735-A-90	17	35	90
G38608	K-OSN-1635-A-90-US	16	35	90
G38681	K-OSN-1730-A-60-US	17	30	60
G38434	K-OSN-1735-A-60-US	17	35	60
G38435	K-OSN-1735-A-90-US	17	35	90
B Bevel				
Heavy Wall				
G30430	K-OSN-1630-B-60	16	30	60
G30447	K-OSN-1630-B-90	16	30	90
G30432	K-OSN-1635-B-60	16	35	60
G29617	K-OSN-1635-B-90	16	35	90
G29618	K-OSN-1730-B-60	17	30	60
G30460	K-OSN-1730-B-90	17	30	90
G29625	K-OSN-1735-B-60	17	35	60
G29626	K-OSN-1735-B-90	17	35	90
G38427	K-OSN-1630-B-60-US	16	30	60

Use reference part numbers with the suffix "-US" if ordering in the U.S., Canada, or Latin America.

Ova-Stiff™ EchoTip® Single Lumen Ovum Aspiration Needle (cont.)

Order Number	Reference Part Number	Needle gage	Needle Length cm	Aspiration Line Length cm
G38428	K-OSN-1630-B-90-US	16	30	90
G38429	K-OSN-1635-B-60-US	16	35	60
G38430	K-OSN-1635-B-90-US	16	35	90
G38433	K-OSN-1730-B-60-US	17	30	60
G48363	K-OSN-1730-B-90-US	17	30	90
G38436	K-OSN-1735-B-60-US	17	35	60
G38437	K-OSN-1735-B-90-US	17	35	90

Use reference part numbers with the suffix "-US" if ordering in the U.S., Canada, or Latin America.


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Single Lumen Ovum Aspiration Needle

Used for ovum aspiration.

- The needle features a B bevel.
- The design of the handle facilitates the rotation of the needle during aspiration. The thumb notch indicates the orientation of the bevel.
- The Luer lock proximal hub accepts syringes.
- EchoTip® technology enhances the visualization of the needle tip under ultrasound.
- Needle only; no aspiration line is included.



Order Number	Reference Part Number	Needle gage	Needle Length cm
B Bevel			
Regular Wall			
G20359	K-OPAA-1635	16	35
G20377	K-OPAA-1730	17	30
G20538	K-OPAA-1735	17	35
G26208	K-OPAA-1835	18	35
G26209	K-OPAA-1935	19	35
G48328	K-OPAA-2030	20	30
G48327	K-OPAA-2035	20	35


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Single Lumen Ovum Aspiration Needle

Used for ovum aspiration.

- The design of the handle facilitates the rotation of the needle during aspiration. The thumb notch indicates the orientation of the bevel.
- EchoTip® technology enhances the visualization of the needle tip under ultrasound.



Order Number	Reference Part Number	Needle gage	Needle Length cm	Aspiration Line Length cm
A Bevel				
Thin Wall				
G55489	K-OPS-1630-T-A-60	16	30	60
G55488	K-OPS-1630-T-A-90	16	30	90
G55249	K-OPS-1635-T-A-60	16	35	60
G55487	K-OPS-1635-T-A-90	16	35	90
G55492	K-OPS-1730-T-A-60	17	30	60
G55252	K-OPS-1735-T-A-90	17	35	90
B Bevel				
Regular Wall				
G55495	K-OPS-1630-R-B-60	16	30	60
G55494	K-OPS-1635-R-B-60	16	35	60
G55493	K-OPS-1730-R-B-60	17	30	60
G55250	K-OPS-1735-R-B-60	17	35	60
G55497	K-OPS-1735-R-B-90	17	35	90
Heavy Wall				
G55483	K-OPS-1630-H-B-60	16	30	60
G55481	K-OPS-1630-H-B-90	16	30	90
G55484	K-OPS-1635-H-B-60	16	35	60
G55248	K-OPS-1635-H-B-90	16	35	90
G55486	K-OPS-1735-H-B-60	17	35	60
G55482	K-OPS-1735-H-B-90	17	35	90


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Single Lumen Ovum Aspiration Needle


Used for ovum aspiration.

- The needle features a B bevel.
- The design of the handle facilitates the rotation of the needle during aspiration. The thumb notch indicates the orientation of the bevel.
- EchoTip® technology enhances the visualization of the needle tip under ultrasound.



Order Number	Reference Part Number	Needle gage	Needle Length cm	Aspiration Line Length cm
B Bevel				
Regular Wall				
G55479	K-OSN-1630-R-B-90	16	30	90
G55478	K-OSN-1635-R-B-90	16	35	90
G55477	K-OSN-1735-R-B-90	17	35	90


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Small Gage Ova-Stiff™ Ovum Aspiration Needle

Used for laparoscopic or ultrasound-guided transvaginal aspiration of oocytes from ovarian follicles.

- The stopper design accepts the tip of a Luer slip syringe, facilitating follicle flushing.
- The EchoTip® enhances the visualization of needle tip when used under ultrasound.



Order Number	Reference Part Number	Needle gage	Needle Length cm	Aspiration Line Length cm
A Bevel				
Regular Wall				
G50735	K-OPS-2130	21	30	60
G50736	K-OPS-2135	21	35	60
Thin Wall				
G19715	K-OPS-6035-RWH-ET	18	35	60
G26182	K-OPS-7030-RWH-ET	19	30	60
G26638	K-OPS-7035-RWH-ET	19	35	60
G26631	K-OPS-8035-RWH-ET	20	35	60
B Bevel				
Thin Wall				
G20943	K-OPS-6035-RWH-B-ET	18	35	60
G26788	K-OPS-7035-RWH-B-ET	19	35	60
G26056	K-OPS-8035-RWH-B-ET	20	35	60

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.