



Intrauterine insemination

We manufacture a range of catheters specifically designed for performing intrauterine insemination (IUI). Our comprehensive IUI catheter line accommodates physician preferences with various configurations, sizes, and rigidities. The catheters have clear depth markers and malleable tips designed to accommodate a broad range of anatomical variations.



Insemi-Form® Catheter

The Insemi-Form Catheter can simplify IUI in a variety of cases where uterine access is restricted by cervical or uterine anomalies. Its pliant cannula provides stiffness to help placement while depth markers provide guidance, and a rounded tip has the potential to reduce trauma to the tissue.

Intrauterine insemination


Insemi-Cath® Insemination Catheter52

Insemi-Form® Catheter52

Shepard Intrauterine Insemination Set53

Soft-Pass™ Coaxial Insemination Catheter54


Soft-Pass™ Insemination Catheter54



Insemi-Cath® Insemination Catheter

Used for the introduction of washed spermatozoa into the uterine cavity.

- The Insemi-Cath Insemination Catheter has a flared proximal end that is designed to affix to a standard Luer slip syringe.
- The catheter is made from nontoxic materials.
- The Insemi-Cath Insemination Catheter is sold in boxes of 20.



Order Number	Reference Part Number	Fr	Length cm
Natural Camber Distal Tip and Open End, Without Positioner			
G16707	J-IUIC-351300	3.5	13
Preset Curved Distal Tip, Closed End, Notch and Adjustable Silicone Positioner			
G17343	J-IUIC-351351	3.5	13
Preset Curved Distal Tip, Open End and Adjustable Silicone Positioner			
G16786	J-IUIC-351304	3.5	13

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Insemi-Form® Catheter

Used for the introduction of washed spermatozoa into the uterine cavity.

- The Insemi-Form Catheter is designed for difficult cases where uterine access is restricted by cervical or uterine anomalies.
- The malleable tip can be formed to any desired curve.
- The pliant cannula adds stiffness, and the catheter features sideport infusion.
- The Insemi-Form Catheter is sold in boxes of 20.



Order Number	Reference Part Number	Outer Sleeve Fr	Outer Sleeve Length cm	Malleable Tip Fr	Malleable Tip Length cm
G18903	J-IFC-090017	9.0	12	6.0	5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Shepard Intrauterine Insemination Set


Used for the introduction of washed spermatozoa into the uterine cavity.

- The adjustable positioner allows the catheter to be placed optimally within the uterus.
- The Shepard Intrauterine Insemination Set is sold in boxes of 20.



Order Number	Reference Part Number	Fr	Length cm
Curved Tip with Stainless Steel Inner Cannula			
G16465	J-IUIE-542009	5.4	20
G16454	J-IUIE-541009	5.4	10
Modified Shepard Intrauterine Insemination Catheters			
Straight Tip without Inner Cannula			
G16468	J-IUIE-542030	5.4	20
Curved Tip without Inner Cannula			
G16469	J-IUIE-542031	5.4	20
G16458	J-IUIE-541031	5.4	10
Curved Tip with Inner Cannula			
G16463	J-IUIE-542001	5.4	20


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Soft-Pass™ Coaxial Insemination Catheter


Used for the introduction of washed spermatozoa into the uterine cavity.

- The soft, flexible inner catheter eases placement and promotes patient comfort.
- The rigid outer catheter provides enhanced catheter control and positioning.
- The coaxial insemination catheter consists of a polyethylene guiding catheter and a polyolefin inner catheter.



Order Number	Reference Part Number	Outer Catheter Fr	Outer Catheter Length cm	Inner Catheter Fr	Inner Catheter Length cm
G18190	J-SPI-068012-CX	6.8	12	4.4	19.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Soft-Pass™ Insemination Catheter

Used for the introduction of washed spermatozoa into the uterine cavity.

- The soft, flexible inner catheter eases placement and promotes patient comfort.
- The rigid outer catheter provides enhanced catheter control and positioning.
- The one-piece insemination catheter has a polyethylene outer sleeve and polyolefin catheter tip.



Order Number	Reference Part Number	Outer Sleeve Fr	Outer Sleeve Length cm	Catheter Tip Fr	Catheter Tip Length cm
G17912	J-SPI-068015	6.8	10	4.4	5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.