

Access the papilla with confidence.



D.A.S.H.[®]

SPHINCTEROTOME WITH DOMETIP[®]

Well known for its naturally contoured DomeTip[®] and reliable orientation, the D.A.S.H. Sphincterotome is designed to aid navigation through the papilla and the fronds, potentially reducing some risks associated with sphincterotomy. And, with the integrated hub, you have the ability to inject while the wire is in the lumen, giving you the efficiency and procedural control you've come to expect from D.A.S.H.



D.A.S.H.®

SPHINCTEROTOME WITH DOMETIP®

Used for cannulation of the ductal system and for sphincterotomy.

- DomeTip designed with anatomy in mind to provide access with potentially reduced trauma
- Custom 3-D forming wire ensures reliable cutting-wire orientation
- Integrated hub provides triple-lumen functionality on a double-lumen platform, giving you the ability to inject while the wire is in the lumen
- Pre-attached racetrack for procedural efficiency and enhanced ergonomics
- Preloaded option for convenience with a .021", .025" or .035" wire guide
- Gold, silver and bronze depth reference markings endoscopically identify relative position of the 25 mm monofilament cutting wire within the duct

Order Number	Reference Part Number	Catheter Size Fr	Cutting Wire Length mm	Wire Guide Diameter inch	Wire Guide Length cm	Color Connection
G26850	DASH-21	5.5	25	.021*	n/a	Black
G25224	DASH-1	6	25	.025**	n/a	Blue
G54193	DASH-35	6.5	25	.035**	n/a	Purple

Preloaded with Tracer Metro® Direct Wire Guide

G25790	DASH-21-480	5.5	25	.021	480	Black
G25225	DASH-260	6	25	.025	260	Blue
G25227	DASH-480	6	25	.025	480	Blue
G25666	DASH-35-480	6.5	25	.035	480	Purple

Preloaded with Acrobat™ Calibrated Tip Wire Guide

G34779	DASH-ACRO-25-450	6	25	.025	450	Blue
G34780	DASH-ACRO-35-260	6.5	25	.035	260	Purple
G34778	DASH-ACRO-35-450	6.5	25	.035	450	Purple

Minimum accessory channel 2.8 mm. *Recommend Tracer Metro® Direct Wire Guide, sold separately. **Recommend Acrobat™ Wire Guide, sold separately. Active cord available separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4